


## Recombinant CD43 Antibody / Rabbit Monoclonal [clone SPN/2049R] (V3651)

Catalog No.	Formulation	Size
V3651-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V3651-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V3651SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug
V3651IHC-7ML	Prediluted in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide; *For IHC use only*	7 ml


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG, kappa
<b>Clone Name</b>	SPN/2049R
<b>Purity</b>	Protein A affinity chromatography
<b>UniProt</b>	P16150
<b>Localization</b>	Cell surface
<b>Applications</b>	Western Blot : 1-2ug/ml Flow Cytometry : 1-2ug/10 <sup>6</sup> cells Immunofluorescence : 1-2ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
<b>Limitations</b>	This recombinant CD43 antibody is available for research use only.




IHC testing of FFPE human tonsil with recombinant CD43 antibody (clone SPN/2049R). Required HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min followed by cooling at RT for 20 min.




SDS-PAGE analysis of purified, BSA-free recombinant CD43 antibody (clone SPN/2049R) as confirmation of integrity and purity.



Immunofluorescence staining of human K562 cells with recombinant CD43 antibody (green, clone SPN/2049R) and NucSpot (red).



Western blot testing of human K562 cell lysate with recombinant CD43 antibody (clone SPN/2049R). Expected molecular weight 45-135 kDa depending on glycosylation level.



FACS staining of PFA-fixed human K562 cells with recombinant CD43 antibody (clone SPN/2049R). Red=isotype control, Blue= recombinant CD43 antibody.

## Description

It recognizes a cell surface glycoprotein of 95/115/135kDa (depending upon the extent of glycosylation), identified as CD43. 70-90% of T-cell lymphomas and from 22-37% of B-cell lymphomas express CD43. No reactivity has been observed with reactive B-cells. So a B-lineage population that co-expresses CD43 is highly likely to be a malignant lymphoma, especially a low-grade lymphoma, rather than a reactive B-cell population. When CD43 antibody is used in combination with anti-CD20, effective immunophenotyping of the lymphomas in formalin-fixed tissues can be obtained. Co-staining of a lymphoid infiltrate with anti-CD20 and anti-CD43 argues against a reactive process and favors a diagnosis of lymphoma.

## Application Notes

Optimal dilution of the recombinant CD43 antibody should be determined by the researcher.

1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

## Immunogen

Recombinant human protein was used as the immunogen for the recombinant CD43 antibody.

## Storage

Store the recombinant CD43 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).