

# Recombinant HLA-DQ Antibody (MHC II) [clone rSPV-L3] (V7843)

Catalog No.	Formulation	Size
V7843-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V7843-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V7843SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

### Recombinant MOUSE MONOCLONAL

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Mouse Monoclonal
Isotype	Mouse IgG1, kappa
Clone Name	rSPV-L3
Purity	Protein G affinity chromatography
UniProt	P01908, P01909, P01920
Localization	Cell surface
Applications	Flow Cytometry : 1-2ug/10^6 cells in 0.1ml Immunofluorescence : 1-2ug/ml
Limitations	This recombinant HLA-DQ antibody is available for research use only.



## **Description**

Recognizes a DQ antigen, which is a dimer of 60kDa. The class II molecule is a heterodimer consisting of an alpha (DQA)

and a beta chain (DQB), both anchored in the membrane. It plays a central role in the immune system by presenting peptides derived from extracellular proteins. Class II molecules are expressed in antigen presenting cells (APC: B Lymphocytes, dendritic cells, macrophages). The alpha chain is approximately 33-35kDa. It is encoded by 5 exons; exon 1 encodes the leader peptide, exons 2 and 3 encode the two extracellular domains, and exon 4 encodes the transmembrane domain and the cytoplasmic tail. Within the DQ molecule both the alpha chain and the beta chain contain the polymorphisms specifying the peptide binding specificities, resulting in up to four different molecules. Typing for these polymorphisms is routinely done for bone marrow transplantation. This MAb strongly blocks cytotoxicity activity of T4-positive cytotoxic T cell clones.

### **Application Notes**

Optimal dilution of the recombinant HLA-DQ antibody should be determined by the researcher.

### **Immunogen**

The T4-positive CTL clone HG-38 was used as the immunogen for the recombinant HLA-DQ antibody.

#### **Storage**

Store the recombinant HLA-DQ antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).