


# Recombinant HSP70 Antibody [clone RM432] (R20447)

Catalog No.	Formulation	Size
R20447-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


## Recombinant RABBIT MONOCLONAL

## **Bulk quote request**


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM432
Purity	Protein A purified from animal origin-free supernatant
UniProt	P0DMV8
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE): 1:100-1:500 Western Blot: 1:2000-1:5000
Limitations	This recombinant HSP70 antibody is available for research use only.



IHC staining of FFPE human breast cancer tissue with recombinant HSP70 antibody at 1:200.



IHC staining of FFPE human liver cancer tissue with recombinant HSP70 antibody at 1:500.



## **Description**

This antibody reacts to human HSP70/Heat shock 70 kDa protein.

### **Application Notes**

The stated application concentrations are suggested starting points. Titration of the recombinant HSP70 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### **Immunogen**

A recombinant human HSP70 protein was used as the immunogen for the recombinant HSP70 antibody.

### **Storage**

Store the recombinant HSP70 antibody at -20oC.