


Recombinant p120-Catenin Antibody / CTNND1 / Catenin delta 1 [clone RM353] (R20375)

Catalog No.	Formulation	Size
R20375-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM353
Purity	Protein A purified from animal origin-free supernatant
UniProt	O60716
Localization	Cell Surface, Cytoplasm
Applications	Immunohistochemistry (FFPE) : 1:1000-1:5000 Western Blot : 1:500-1:1000
Limitations	This recombinant p120-Catenin antibody is available for research use only.



Western blot testing of human A431 cell lysate with recombinant p120-Catenin antibody at 1:500 dilution. Predicted molecular weight of ~108 kDa but commonly observed at 66-120 kDa (multiple isoforms).



IHC staining of FFPE human kidney tissue with recombinant p120-Catenin antibody at 1:5000.

Description

This antibody reacts human p120 catenin (Catenin delta-1).

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant p120-Catenin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of human p120 catenin (Catenin delta-1) was used as the immunogen for the recombinant p120-Catenin antibody.

Storage

Store the recombinant p120-Catenin antibody at -20oC.