


## Recombinant Progesterone Antibody [clone RM434] (R20449)

Catalog No.	Formulation	Size
R20449-0.1ML	1mg/ml in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM434
Purity	Protein A purified from animal origin-free supernatant
Applications	Competitive ELISA : 1ng/well-10ng/well
Limitations	This recombinant Progesterone antibody is available for research use only.



Representative standard curve for Progesterone Competitive ELISA.

### Description

This antibody reacts to Progesterone.

### Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Progesterone antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

KLH-conjugated Progesterone was used as the immunogen for the recombinant Progesterone antibody.

## Storage

Store the recombinant Progesterone antibody at -20oC.