


## Recombinant Prolactin Antibody / PRL [clone PRL/4906R] (V8869)

Catalog No.	Formulation	Size
V8869-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V8869-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V8869SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	PRL/4906R
<b>Purity</b>	Protein A/G affinity
<b>UniProt</b>	IDP01236
<b>Localization</b>	Secreted
<b>Applications</b>	Immunohistochemistry (FFPE) : 1-2ug/ml
<b>Limitations</b>	This recombinant Prolactin antibody is available for research use only.



IHC staining of FFPE human pituitary tissue with recombinant Prolactin antibody (clone PRL/4906R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE human pituitary tissue with recombinant Prolactin antibody (clone PRL/4906R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



Negative control: IHC staining of FFPE human brain tissue using recombinant Prolactin antibody (clone PRL/4906R) at 2ug/ml in PBS for 30min RT. HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.

## Description

Prolactin is a growth factor that is secreted by the anterior pituitary. It is necessary for the proliferation and differentiation of the mammary glands. Prolactin is useful in the classification of pituitary tumors and study of pituitary disease. It also plays a role in the development of mammary cancer, functioning dually as a mitogen and a differentiating agent.

## Application Notes

Optimal dilution of the recombinant Prolactin antibody should be determined by the researcher.

## Immunogen

A portion of amino acids 63-201 was used as the immunogen for the recombinant Prolactin antibody.

## Storage

Aliquot the recombinant Prolactin antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.