


## Recombinant SOX2 Antibody [clone SOX2/3811R] (V7864)

Catalog No.	Formulation	Size
V7864-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V7864-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V7864SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	SOX2/3811R
Purity	Protein A affinity chromatography
UniProt	P48431
Localization	Nuclear
Applications	ELISA (order BSA-free Format For Coating) : Western Blot : 1-2ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This recombinant SOX2 antibody is available for research use only.



IHC staining of FFPE human cervical carcinoma with recombinant SOX2 antibody (clone SOX2/3811R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min and allow to cool before testing.



IHC staining of FFPE human cervical carcinoma with recombinant SOX2 antibody (clone SOX2/3811R). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 10-20 min and allow to cool before testing.

## Description

SOX2 is required for stem cell maintenance in the central nervous system, and it also regulates gene expression in the stomach. SOX2 is necessary for regulating multiple transcription factors that affect Oct 3/4 expression. An essential function of SOX2 is to stabilize embryonic stem cells in a pluripotent state by maintaining the requisite level of Oct 3/4 expression. Reportedly, SOX2 is associated with aggressive phenotypes of breast, head and neck, gastric, colorectal, bladder, and small cell lung cancers. However, SOX2 is expressed in a high percentage of lung squamous cell carcinomas and has been shown to be an independent favorable prognostic marker.

## Application Notes

Optimal dilution of the recombinant SOX2 antibody should be determined by the researcher.

## Immunogen

A recombinant human partial protein (amino acids 176-305) was used as the immunogen for the recombinant SOX2 antibody.

## Storage

Store the recombinant SOX2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).