


Recombinant Testosterone Antibody [clone RM435] (R20450)

Catalog No.	Formulation	Size
R20450-0.1ML	1mg/ml in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM435
Purity	Protein A purified from animal origin-free supernatant
Applications	Competitive ELISA : 1ng/well-10ng/well
Limitations	This recombinant Testosterone antibody is available for research use only.



Representative standard curve for Testosterone Competitive ELISA.

Description

This antibody reacts to Testosterone.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Testosterone antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

KLH-conjugated Testosterone was used as the immunogen for the recombinant Testosterone antibody.

Storage

Store the recombinant Testosterone antibody at -20oC.