


Recombinant TLR2 Antibody [clone TLR2/3894R] (V8393)

Catalog No.	Formulation	Size
V8393-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V8393-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V8393SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	TLR2/3894R
Purity	Protein A affinity chromatography
UniProt	O60603
Localization	Cell surface
Applications	ELISA : For coating order antibody without BSA
Limitations	This recombinant TLR2 antibody is available for research use only.



SDS-PAGE analysis of purified, BSA-free recombinant TLR2 antibody (clone TLR2/3894R) as confirmation of integrity and purity.

Description

This MAb reacts with human Toll-like receptor 2 (TLR2). It is a member of the Toll-like receptor (TLR) family, which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from *Drosophila* to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. This gene is expressed most abundantly in peripheral blood leukocytes, and mediates host response to Gram-positive bacteria and yeast via stimulation of NF-kappaB.

Application Notes

Optimal dilution of the recombinant TLR2 antibody should be determined by the researcher.

Immunogen

Recombinant full-length human TLR2 protein was used as the immunogen for the recombinant TLR2 antibody.

Storage

Store the recombinant TLR2 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).