


RGS13 Antibody (R33428)

Catalog No.	Formulation	Size
R33428-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6003
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RGS13 antibody is available for research use only.



Western blot testing of peripheral blood mononucleocyte lysate with RGS13 antibody at 0.5ug/ml. Predicted molecular weight: ~19kDa.

Description

Additional name(s) for this target protein: Regulator of G-protein signaling 13

Application Notes

Optimal dilution of the RGS13 antibody should be determined by the researcher.

Immunogen

Amino acids DESKRPPSNLTLEE were used as the immunogen for this RGS13 antibody.

Storage

Aliquot and store the RGS13 antibody at -20°C.