


RGS18 Antibody (R33286)

Catalog No.	Formulation	Size
R33286-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64407
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RGS18 antibody is available for research use only.



Western blot testing of human peripheral blood mononucleocyte lysate with RGS18 antibody at 1ug/ml. Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Regulator of G-protein signalling 18

Application Notes

Optimal dilution of the RGS18 antibody should be determined by the researcher.

Immunogen

Amino acids CNEFQDVQSDVAIWL were used as the immunogen for this RGS18 antibody.

Storage

Aliquot and store the RGS18 antibody at -20oC.