


## RHBG Antibody (R35066)

Catalog No.	Formulation	Size
R35066-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	58176
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This RHBG antibody is available for research use only.



Western blot testing of NSO lysate with RHBG antibody at 1ug/ml. Predicted molecular weight: ~50/55kDa (unmodified/glycosylated).

## Description

Additional name(s) for this target protein: Rh family B glycoprotein

## Application Notes

Optimal dilution of the RHBG antibody should be determined by the researcher.

## Immunogen

Amino acids PLRGGESDTR were used as the immunogen for this RHBG antibody.

## Storage

Aliquot and store the RHBG antibody at -20oC.

## References (1)