


RHCG Antibody (R33974)

Catalog No.	Formulation	Size
R33974-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51458
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RHCG antibody is available for research use only.



Western blot testing of human testis lysate with RHCG antibody at 1ug/ml. Predicted molecular weight: ~53kDa, routinely observed at ~60kDa.

Description

Additional name(s) for this target protein: Rh family, C glycoprotein

Application Notes

Optimal dilution of the RHCG antibody should be determined by the researcher.

Immunogen

Amino acids ERQNSVYQSD were used as the immunogen for this RHCG antibody.

Storage

Aliquot and store the RHCG antibody at -20oC.