


RHOGDI Antibody / ARHGDIA (R35344)

Catalog No.	Formulation	Size
R35344-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	396
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This RHOGDI antibody is available for research use only.



Western blot testing of mouse lung lysate with RHOGDI antibody at 1ug/ml. Predicted molecular weight: ~23 kDa, routinely observed at ~29 kDa.

Description

Additional name(s) for this target protein: Rho GDP dissociation inhibitor alpha; ARHGDIA

Application Notes

Optimal dilution of the RHOGDI antibody should be determined by the researcher.

Immunogen

Amino acids QLAQIAAENE were used as the immunogen for this RHOGDI antibody.

Storage

Aliquot and store the RHOGDI antibody at -20oC.