


RIC8A Antibody (R33280)

Catalog No.	Formulation	Size
R33280-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	60626
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RIC8A antibody is available for research use only.



Western blot testing of human brain lysate with RIC8A antibody at 1ug/ml. Predicted molecular weight: ~60kDa, observed here at ~70kDa.

Description

Additional name(s) for this target protein: Resistance to inhibitors of cholinesterase 8a

Application Notes

Optimal dilution of the RIC8A antibody should be determined by the researcher.

Immunogen

Amino acids CLSSDPDSDPD were used as the immunogen for this RIC8A antibody.

Storage

Aliquot and store the RIC8A antibody at -20oC.