


## RICTOR Antibody (R34597)

Catalog No.	Formulation	Size
R34597-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	253260
Applications	Immunohistochemistry (FFPE) : 2-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This RICTOR antibody is available for research use only.



IHC testing of FFPE human kidney with RICTOR antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: Rapamycin-insensitive companion of mTOR

### Application Notes

Optimal dilution of the RICTOR antibody should be determined by the researcher.

## **Immunogen**

Amino acids KQPIVD TSAES were used as the immunogen for this RICTOR antibody.

## **Storage**

Aliquot and store the RICTOR antibody at -20oC.