


RKIP Antibody / PBP / PEBP1 Biotin Conjugate (R33567BTN)

Catalog No.	Formulation	Size
R33567BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5037
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This RKIP antibody is available for research use only.



Western blot testing of human cerebellum lysate with biotinylated RKIP antibody at 0.1ug/ml. Predicted molecular weight: ~21kDa.

Description

Additional name(s) for this target protein: RAF kinase inhibitor protein; PEBP1

Application Notes

Optimal dilution of the RKIP antibody should be determined by the researcher.

Immunogen

Amino acids DPDAPSRKDPKYRE were used as the immunogen for this RKIP antibody.

Storage

Aliquot and store the RKIP antibody at -20oC.