


RKIP Antibody / PBP / PEBP1 (R33567)


Catalog No.	Formulation	Size
R33567-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5037
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This RKIP antibody is available for research use only.



Western blot testing of human prostate lysate with RKIP antibody at 0.01ug/ml.
Predicted molecular weight: ~21kDa.



IHC staining of FFPE human adrenal gland with RKIP antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: RAF kinase inhibitor protein; PEBP1

Application Notes

Optimal dilution of the RKIP antibody should be determined by the researcher.

Immunogen

Amino acids DPDAPSRKDPKYRE were used as the immunogen for this RKIP antibody.

Storage

Aliquot and store the RKIP antibody at -20oC.