


RNF115 Antibody (R33588)

Catalog No.	Formulation	Size
R33588-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27246
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RNF115 antibody is available for research use only.



Western blot testing of HEPG2 cell lysate with RNF115 antibody at 0.3ug/ml. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: Ring finger protein 115; breast cancer-associated gene 2

Application Notes

Optimal dilution of the RNF115 antibody should be determined by the researcher.

Immunogen

Amino acids DQDNRANERGHQTH were used as the immunogen for this RNF115 antibody.

Storage

Aliquot and store the RNF115 antibody at -20oC.