


RNF13 Antibody (R34215)

Catalog No.	Formulation	Size
R34215-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11342
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This RNF13 antibody is available for research use only.



Description

Additional name(s) for this target protein: Ring finger protein 13

Application Notes

Optimal dilution of the RNF13 antibody should be determined by the researcher.

Immunogen

Amino acids HKFKKGDEYDVC were used as the immunogen for this RNF13 antibody. **Storage** Aliquot and store the RNF13 antibody at -20oC.