

RNF13 Antibody (R35077)

Catalog No.	Formulation	Size
R35077-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11342
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml (2) ELISA (peptide) LOD : 1:1000
Limitations	This RNF13 antibody is available for research use only.



Description

Additional name(s) for this target protein: Ring finger protein 13, RZF

Application Notes

Optimal dilution of the RNF13 antibody should be determined by the researcher.

1. This RNF13 antibody will detect isoforms 1 and 3.
2. This antibody stained islets but not in ducts and acini of Human normal pancreas. However, in pancreatic ductal adenocarcinoma staining was in the ducts and acini.

Immunogen

Amino acids PNGERDYNIAITV were used as the immunogen for this RNF13 antibody.

Storage

Aliquot and store the RNF13 antibody at -20°C.