


RNF213 Antibody (R34659)

Catalog No.	Formulation	Size
R34659-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57674
Applications	Immunohistochemistry (FFPE) : 2.5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This RNF213 antibody is available for research use only.



IHC testing of FFPE human kidney with RNF213 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining. Rare nucleoli staining is seen.

Description

Additional name(s) for this target protein: Ring finger protein 213, C17orf27

Application Notes

Optimal dilution of the RNF213 antibody should be determined by the researcher.

Immunogen

Amino acids KTA AVLKWNRE were used as the immunogen for this RNF213 antibody.

Storage

Aliquot and store the RNF213 antibody at -20oC.