


RNF31 Antibody (R33762)

Catalog No.	Formulation	Size
R33762-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55072
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RNF31 antibody is available for research use only.



Western blot testing of HepG2 lysate with RNF31 antibody at 2ug/ml. Predicted molecular weight: ~120kDa.

Description

Additional name(s) for this target protein: Ring finger protein 31; ZIBRA

Application Notes

Optimal dilution of the RNF31 antibody should be determined by the researcher.

Immunogen

Amino acids EEVPLGQSIPRRRK were used as the immunogen for this RNF31 antibody.

Storage

Aliquot and store the RNF31 antibody at -20oC.