


RNF7 Antibody (R33370)

Catalog No.	Formulation	Size
R33370-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9616
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This RNF7 antibody is available for research use only.



Western blot testing of human heart lysate with RNF7 antibody at 1ug/ml. Predicted molecular weight: ~13kDa.

Description

Additional name(s) for this target protein: Ring finger protein 7

Application Notes

Optimal dilution of the RNF7 antibody should be determined by the researcher.

1. This RNF7 antibody is specific for isoform 1.

Immunogen

Amino acids DACLRCQAENKQE were used as the immunogen for this RNF7 antibody.

Storage

Aliquot and store the RNF7 antibody at -20oC.