


## RNF7 Antibody (R33370)

Catalog No.	Formulation	Size
R33370-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9616
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This RNF7 antibody is available for research use only.



Western blot testing of human heart lysate with RNF7 antibody at 1ug/ml. Predicted molecular weight: ~13kDa.

## Description

Additional name(s) for this target protein: Ring finger protein 7

## Application Notes

Optimal dilution of the RNF7 antibody should be determined by the researcher.

1. This RNF7 antibody is specific for isoform 1.

## **Immunogen**

Amino acids DACLRCQAENKQE were used as the immunogen for this RNF7 antibody.

## **Storage**

Aliquot and store the RNF7 antibody at -20oC.