


## RNF8 Antibody (R34083)

Catalog No.	Formulation	Size
R34083-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	9025
<b>Applications</b>	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:64000
<b>Limitations</b>	This RNF8 antibody is available for research use only.



Western blot testing of human lung lysate with RNF8 antibody at 0.1ug/ml. Predicted molecular weight ~56 kDa but may be observed at larger sizes due to ubiquitination.

### Description

Additional name(s) for this target protein: Ring finger protein 8

### Application Notes

Optimal dilution of the RNF8 antibody should be determined by the researcher.

### Immunogen

Amino acids GEPGFFVTGDRAG were used as the immunogen for this RNF8 antibody.

## **Storage**

Aliquot and store the RNF8 antibody at -20°C.