


## RNF8 Antibody (R34083)

Catalog No.	Formulation	Size
R34083-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9025
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RNF8 antibody is available for research use only.



Western blot testing of human lung lysate with RNF8 antibody at 0.1ug/ml. Predicted molecular weight ~56 kDa but may be observed at larger sizes due to ubiquitination.

## Description

Additional name(s) for this target protein: Ring finger protein 8

## Application Notes

Optimal dilution of the RNF8 antibody should be determined by the researcher.

## Immunogen

Amino acids GEPGFFVTGDRAQ were used as the immunogen for this RNF8 antibody.

## Storage

Aliquot and store the RNF8 antibody at -20oC.