


RNF8 Antibody (R34083)

Catalog No.	Formulation	Size
R34083-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9025
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RNF8 antibody is available for research use only.



Western blot testing of human lung lysate with RNF8 antibody at 0.1ug/ml. Predicted molecular weight ~56 kDa but may be observed at larger sizes due to ubiquitination.

Description

Additional name(s) for this target protein: Ring finger protein 8

Application Notes

Optimal dilution of the RNF8 antibody should be determined by the researcher.

Immunogen

Amino acids GEPGFFVTGDRAG were used as the immunogen for this RNF8 antibody.

Storage

Aliquot and store the RNF8 antibody at -20oC.