


RORC Antibody (R33246)

Catalog No.	Formulation	Size
R33246-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6097
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RORC antibody is available for research use only.



Western blot testing of human colon lysate with RORC antibody at 0.1ug/ml. Predicted molecular weight: ~56 kDa.

Description

Additional name(s) for this target protein: RAR-related orphan receptor C

Application Notes

Optimal dilution of the RORC antibody should be determined by the researcher.

Immunogen

Amino acids CHLEYSERGKAE were used as the immunogen for this RORC antibody.

Storage

Aliquot and store the RORC antibody at -20oC.