


# **RORC Antibody (R33246)**

Catalog No.	Formulation	Size
R33246-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6097
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RORC antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: RAR-related orphan receptor C

## **Application Notes**

Optimal dilution of the RORC antibody should be determined by the researcher.

#### **Immunogen**

Amino acids CHLEYSPERGKAE were used as the immunogen for this RORC antibody. **Storage** Aliquot and store the RORC antibody at -20oC.