


RPA70 Antibody (R33266)


Catalog No.	Formulation	Size
R33266-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6117
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This RPA70 antibody is available for research use only.



Western blot testing of HeLa lysate with RPA70 antibody at 0.3ug/ml. Expected molecular weight: ~70kDa. The ~80kDa band is likely acetylated RPA70.



Western blot testing of NIH3T3 (A), 3T3-L1 (B), mouse lymph node (C) and NSO (D) lysates with RPA70 antibody at 0.5ug/ml. Expected molecular weight: ~70kDa.

Description

Additional name(s) for this target protein: Replication protein A1, 70kDa; RPA1

Application Notes

Optimal dilution of the RPA70 antibody should be determined by the researcher.

Immunogen

Amino acids CKSYEDATKITVRSN were used as the immunogen for this RPA70 antibody.

Storage

Aliquot and store the RPA70 antibody at -20oC.