

## RPGRIP1L Antibody / Protein phantom (R33497)

Catalog No.	Formulation	Size
R33497-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23322
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RPGRIP1L antibody is available for research use only.



Western blot testing of human cerebellum lysate with RPGRIP1L antibody at 0.3ug/ml.  
Predicted molecular weight: ~151kDa.

### Description

Additional name(s) for this target protein: RPGRIP1-like, Phantom, FTM

### Application Notes

Optimal dilution of the RPGRIP1L antibody should be determined by the researcher.

## Immunogen

Amino acids DKENNKAKRDILK were used as the immunogen for this RPGRIP1L antibody.

## Storage

Aliquot and store the RPGRIP1L antibody at -20oC.