


# RPL11 Antibody / 60S ribosomal protein L11 (RQ7875)


Catalog No.	Formulation	Size
RQ7875	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

## **Bulk quote request**


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P62913
Localization	Cytoplasm, nucleus, nucleolus
Applications	Western Blot: 0.5-1ug/ml Immunohistochemistry (FFPE): 2-5ug/ml Immunofluorescence (FFPE): 5ug/ml Flow Cytometry: 1-3ug/million cells Direct ELISA: 0.1-0.5ug/ml
Limitations	This RPL11 antibody is available for research use only.




IHC staining of FFPE mouse brain tissue with RPL11 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE rat brain tissue with RPL11 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE placental tissue with RPL11 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Immunofluorescent staining of FFPE human MCF7 cells with RPL11 antibody (red) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of 1) human A549, 2) human Jurkat, 3) human Raji, 4) human HeLa, 5) mouse liver and 6) mouse NIH 3T3 cell lysate with RPL11 antibody. Predicted molecular weight ~18 kDa.



Flow cytometry testing of human ThP-1 cells with RPL11 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= RPL11 antibody.

### **Description**

60S ribosomal protein L11 is a protein that in humans is encoded by the RPL11 gene. Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L5P family of ribosomal proteins. It is located in the cytoplasm. The protein probably associates with the 5S rRNA. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

## **Application Notes**

Optimal dilution of the RPL11 antibody should be determined by the researcher.

#### **Immunogen**

E. coli-derived recombinant human protein (amino acids E7-K178) was used as the immunogen for the RPL11 antibody.

#### **Storage**

After reconstitution, the RPL11 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.