


RRAD Antibody (R34275)

Catalog No.	Formulation	Size
R34275-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6236
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RRAD antibody is available for research use only.



Description

Additional name(s) for this target protein: Ras-related associated with diabetes

Application Notes

Optimal dilution of the RRAD antibody should be determined by the researcher.

Immunogen

Amino acids HRRSMPVDERDLQ were used as the immunogen for this RRAD antibody.

Storage

Aliquot and store the RRAD antibody at -20oC.