

RRAD Antibody (R34275)

Catalog No.	Formulation	Size
R34275-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6236
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RRAD antibody is available for research use only.




Figure A): Western blot of HEK293 lysates overexpressing several HA-tagged mouse GTPases including RRAD and compared with an HA-specific antibody. Figure B): Western blot of WT and KO lysates of mouse heart (100ug protein). RRAD antibody was used at 0.5ug/ml.

Description

Additional name(s) for this target protein: Ras-related associated with diabetes

Application Notes

Optimal dilution of the RRAD antibody should be determined by the researcher.

Immunogen

Amino acids HRRSMPVDERDLQ were used as the immunogen for this RRAD antibody. **Storage** Aliquot and store the RRAD antibody at -20oC.