


RSL1D1 Antibody (R35167)

Catalog No.	Formulation	Size
R35167-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26156
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This RSL1D1 antibody is available for research use only.



Western blot testing of K562 lysate with RSL1D1 antibody at 1ug/ml. Predicted molecular weight: ~55kDa, observed here at ~75kDa.

Description

Additional name(s) for this target protein: Ribosomal L1 domain containing 1, CSIG, cellular senescence inhibited gene protein

Application Notes

Optimal dilution of the RSL1D1 antibody should be determined by the researcher.

Immunogen

Amino acids PTARKQLDKEQVRK were used as the immunogen for this RSL1D1 antibody.

Storage

Aliquot and store the RSL1D1 antibody at -20oC.