



## Runx1 Antibody (F53249)

Catalog No.	Formulation	Size
F53249-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53249-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q03347
<b>Applications</b>	Western Blot : 1:2000
<b>Limitations</b>	This Runx1 antibody is available for research use only.



## Description

Runx1 is essential for the development of normal hematopoiesis.

## **Application Notes**

Titration of the Runx1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

A portion of amino acids 300-350 from the mouse protein was used as the immunogen for this Runx1 antibody.

## **Storage**

Aliquot the Runx1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.