

Runx1 Antibody (F53249)


Catalog No.	Formulation	Size
F53249-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F53249-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q03347
Applications	Western Blot : 1:2000
Limitations	This Runx1 antibody is available for research use only.



Western blot testing of Runx1 antibody at 1:2000 dilution + SW620 lysate; Predicted molecular weight: 49 kDa.



Description

Runx1 is essential for the development of normal hematopoiesis.

Application Notes

Titration of the Runx1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 300-350 from the mouse protein was used as the immunogen for this Runx1 antibody.

Storage

Aliquot the Runx1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.