


RXR alpha Antibody (R33316)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33316-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Mouse, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 6256 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This RXR alpha antibody is available for research use only. |



Western blot testing of nuclear HeLa lysate with RXR alpha antibody at 0.3ug/ml.
Predicted molecular weight: ~51kDa but routinely observed at 54~60kDa.

Description

Additional name(s) for this target protein: Retinoid X receptor alpha; RXRA

Application Notes

Optimal dilution of the RXR alpha antibody should be determined by the researcher.

1. This RXR alpha antibody does not react with RXRb or RXRg.

Immunogen

Amino acids CQVNSSLTSPTRGSM were used as the immunogen for this RXR alpha antibody.

Storage

Aliquot and store the RXR alpha antibody at -20oC.

References (1)