


RXRG Antibody (R35739)

Catalog No.	Formulation	Size
R35739-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6258
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This RXRG antibody is available for research use only.



Western blot testing of human brain lysate with RXRG antibody at 1ug/ml. Predicted molecular weight: ~51 kDa.

Description

Additional name(s) for this target protein: Retinoid X receptor gamma, RXR gamma

Application Notes

Optimal dilution of the RXRG antibody should be determined by the researcher.

Immunogen

Amino acids RQSRERAEESEAEK were used as the immunogen for this RXRG antibody. No cross reaction with RXRA or RXRB expected.

Storage

Aliquot and store the RXRG antibody at -20oC.