


## S100A8 Antibody (R34437)

Catalog No.	Formulation	Size
R34437-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	6279
<b>Applications</b>	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This S100A8 antibody is available for research use only.



Western blot testing of human lung lysate with S100A8 antibody at 0.03ug/ml. Predicted molecular weight ~11kDa.

### Description

Additional name(s) for this target protein: S100 calcium binding protein A8

### Application Notes

Optimal dilution of the S100A8 antibody should be determined by the researcher.

### Immunogen

Amino acids KGADVWFKELDINT were used as the immunogen for this S100A8 antibody.

## **Storage**

Aliquot and store the S100A8 antibody at -20oC.

## **Alternate Names**

Calprotectin L1L subunit , MRP8, Calgranulin A