


\$100A8 Antibody (R34437)

Catalog No.	Formulation	Size
R34437-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6279
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This S100A8 antibody is available for research use only.



Description

Additional name(s) for this target protein: S100 calcium binding protein A8

Application Notes

Optimal dilution of the S100A8 antibody should be determined by the researcher.

Immunogen

Amino acids KGADVWFKELDINT were used as the immunogen for this S100A8 antibody.

Storage

Aliquot and store the S100A8 antibody at -20oC.

Alternate Names

Calprotectin L1L subunit , MRP8, Calgranulin A