


Sall4 Antibody (R33603)

Catalog No.	Formulation	Size
R33603-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	99377
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Sall4 antibody is available for research use only.



Western blot testing of rat testis lysate with Sall4 antibody at 0.1ug/ml. Predicted molecular weight ~110/68kDa (Isoforms a/b). All three bands are blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Sal-like 4

Application Notes

Optimal dilution of the Sall4 antibody should be determined by the researcher.

1. This Sall4 antibody is specific for isoforms a and b of the protein.

Immunogen

Amino acids DQSKKGKGKPQN were used as the immunogen for this Sall4 antibody.

Storage

Aliquot and store the Sall4 antibody at -20oC.