


SAP97 Antibody (R32206)


Catalog No.	Formulation	Size
R32206	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
UniProt	Q12959
Localization	Cytoplasmic, membrane
Applications	Western Blot : 0.1-0.5ug/ml IHC (FFPE) : 0.5-1ug/ml
Limitations	This SAP97 antibody is available for research use only.



Western blot testing of 1) rat lung, 2) mouse lung, 3) HeLa, 4) MM231, 5) COLO320, 6) A549 and 7) mouse NIH3T3 lysate with SAP97 antibody. Predicted molecular weight: ~103 kDa, routinely observed at ~140 kDa (Ref 1).



IHC testing of FFPE human prostate cancer tissue with SAP97 antibody. HIER: Boil the paraffin sections in pH 6, 10mM citrate buffer for 20 minutes and allow to cool prior to staining.

Description

Disks large homolog 1 (DLG1), also known as synapse-associated protein 97 or SAP97, is a protein that in humans is encoded by the SAP97 gene. SAP97 is expressed throughout the body in epithelial cells, including the kidney and brain. There is some evidence that SAP97 regulates cell-to-cell adhesion during cell death, and may interact with HPV. In the brain, SAP97's function is involved in the trafficking of transmembrane receptors from the ER to the plasma membrane. SAP97's function has been investigated by reducing its expression by knockout or increasing its expression heterologously. Mice in which the SAP97 gene has been knocked out die perinatally, have a cleft palate, and deficiencies in renal function. Overexpression of SAP97 in mammalian neurons leads to increased synaptic strength.

Application Notes

Optimal dilution of the SAP97 antibody should be determined by the researcher.

Immunogen

Amino acids 1-165 of human SAP97 were used as the immunogen for the SAP97 antibody.

Storage

After reconstitution, the SAP97 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.

References (1)