

SAR1B Antibody (R33465)


Catalog No.	Formulation	Size
R33465-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51128
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This SAR1B antibody is available for research use only.



Western blot testing of mouse liver lysate with SAR1B antibody at 0.03ug/ml. Predicted molecular weight: ~22kDa.



Western blot of HEK293 lysate overexpressing SAR1B probed with SAR1B antibody (mock transfection in lane 1). Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: SAR1a gene homolog 2

Application Notes

Optimal dilution of the SAR1B antibody should be determined by the researcher.

Immunogen

Amino acids DHERLLESKEELD were used as the immunogen for this SAR1B antibody.

Storage

Aliquot and store the SAR1B antibody at -20oC.