


## SART3 Antibody / Tip110 (R34754)

Catalog No.	Formulation	Size
R34754-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9733
Applications	Western Blot : 0.01-0.01ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This SART3 antibody is available for research use only.



Western blot testing of Jurkat cell lysate with SART3 antibody at 0.01ug/ml. Predicted molecular weight: ~111kDa, observed here at ~170kDa.

## Description

Additional name(s) for this target protein: Squamous cell carcinoma antigen recognised by T cells 3; Tip110; Tat-interacting protein of 110 kDa; p110 nuclear RNA-binding protein

## Application Notes

Optimal dilution of the SART3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids MSNADFAKLFLRK were used as the immunogen for this SART3 antibody.

## **Storage**

Aliquot and store the SART3 antibody at -20oC.