


SCARB2 Antibody Biotin Conjugate (R34831BTN)

Catalog No.	Formulation	Size
R34831BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	950
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:4000
Limitations	This SCARB2 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with biotinylated SCARB2 antibody at 1ug/ml. Predicted molecular weight: ~54/85kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Scavenger receptor class B, member 2; SR-BII; LIMP2; LIMPII; LGP85; 85 kDa lysosomal sialoglycoprotein scavenger receptor class B, member 2

Application Notes

Optimal dilution of the SCARB2 antibody should be determined by the researcher.

Immunogen

Amino acids NKANIQFGDNGTTIS were used as the immunogen for this SCARB2 antibody.

Storage

Aliquot and store the SCARB2 antibody at -20oC.

References (1)