


## SCG10 Antibody (R33144)

Catalog No.	Formulation	Size
R33144-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Mouse
<b>Predicted Reactivity</b>	Cow, Dog, Human, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	11075
<b>Localization</b>	Cytoplasmic
<b>Applications</b>	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This SCG10 antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with SCG10 antibody at 0.5ug/ml.  
Predicted molecular weight: ~21 kDa.

### Description

Additional name(s) for this target protein: Stathmin-like 2; STMN2

### Application Notes

Optimal dilution of the SCG10 antibody should be determined by the researcher.

## **Immunogen**

Amino acids AKTAMAYKEK were used as the immunogen for this SCG10 antibody.

## **Storage**

Aliquot and store the SCG10 antibody at -20oC.