


SCG10 Antibody (R33144)

Catalog No.	Formulation	Size
R33144-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Dog, Human, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11075
Localization	Cytoplasmic
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SCG10 antibody is available for research use only.



Western blot testing of mouse fetal brain lysate with SCG10 antibody at 0.5ug/ml.
Predicted molecular weight: ~21 kDa.

Description

Additional name(s) for this target protein: Stathmin-like 2; STMN2

Application Notes

Optimal dilution of the SCG10 antibody should be determined by the researcher.

Immunogen

Amino acids AKTAMAYKEK were used as the immunogen for this SCG10 antibody.

Storage

Aliquot and store the SCG10 antibody at -20oC.