


## SCGN Antibody [clone AFGD-19] (RQ4793)

Catalog No.	Formulation	Size
RQ4793	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

<b>Availability</b>	1-2 weeks
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	AFGD-19
<b>Purity</b>	Affinity purified
<b>UniProt</b>	O76038
<b>Applications</b>	Western Blot : 1:500-1:2000
<b>Limitations</b>	This SCGN antibody is available for research use only.



Western blot testing of human fetal brain lysate with SCGN antibody. Predicted molecular weight ~32 kDa.

### Description

Secretagoin is expressed at high levels in the pancreatic islets of Langerhans and to a much lesser extent in the gastrointestinal tract (stomach, small intestine and colon), the adrenal medulla and cortex and the thyroid C-cells.

### Application Notes

Optimal dilution of the SCGN antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide specific to human SCGN was used as the immunogen for the SCGN antibody.

## **Storage**

Store the SCGN antibody at -20oC.